

COVID-19 DETECTION SYSTEM

THEME: DEEP LEARNING, COMPUTER VISION AND
IMAGE PROCESSING

DEVELOPED BY,
SURYA TEJA CHALLAPALLI (9.6 YEARS OF EXP)

SOLUTION LINKS & TECHNOLOGIES

- <https://awschalst1988.notebook.us-east-1.sagemaker.aws/lab>
- <https://github.com/chalst1988/Competitions>





- Amazon Web Services (AWS)
- **Amazon Sage Maker Studio**
- Jupyter Lab
- Notebook
- Python
- Deep Learning, Tensorflow, Keras

ABSTRACT & APPROACH

X-ray images are given in jpeg, jpg and png file formats. Folder structure contains Train and Test folders. Each folder is segregated into Covid and Normal folders as shown below. Train contains 224 images and Test contains 60 images.

- Initially a Convolutional Neural Network was build to train the images with pooling, dropout, flatten and dense layers.
- Secondly results were predicted as 0 and 1 for the given train images. {'Covid':0, 'Normal':1}
- Finally labelling on image were done as “COVID” or “NORMAL” based on the predicted classes.





```
print(model.summary())
```


Model: "sequential_1"

Layer (type)	Output Shape	Param #
=====		
conv2d_1 (Conv2D)	(None, 222, 222, 32)	896
conv2d_2 (Conv2D)	(None, 220, 220, 128)	36992
max_pooling2d_1 (MaxPooling2D)	(None, 110, 110, 128)	0
dropout_1 (Dropout)	(None, 110, 110, 128)	0
conv2d_3 (Conv2D)	(None, 108, 108, 64)	73792
max_pooling2d_2 (MaxPooling2D)	(None, 54, 54, 64)	0
dropout_2 (Dropout)	(None, 54, 54, 64)	0
conv2d_4 (Conv2D)	(None, 52, 52, 128)	73856
max_pooling2d_3 (MaxPooling2D)	(None, 26, 26, 128)	0
dropout_3 (Dropout)	(None, 26, 26, 128)	0
flatten_1 (Flatten)	(None, 86528)	0
dense_1 (Dense)	(None, 64)	5537856
dropout_4 (Dropout)	(None, 64)	0
dense_2 (Dense)	(None, 1)	65
=====		
Total params: 5,723,457		
Trainable params: 5,723,457		
Non-trainable params: 0		
None		


MODEL SUMMARY & HISTORY PLOTS

VISUALIZE TEST X-RAY IMAGES

Prediction of the image: COVID-19



Prediction of the image: NORMAL



Prediction of the image: COVID-19

